Comparison of radial pressure wave therapy and focused shock wave therapy


In recent years, shock wave therapy has established itself as a proven and effective method for treating pain and various musculoskeletal disorders. There are two main forms: radial pressure wave therapy (rESWT) and focused shock wave therapy (fESWT). Although both methods work with mechanical waves, there are fundamental differences in the type of waves generated, their physical properties and the areas of application.

Radial

Radial pressure wave therapy generates its waves using a ballistic principle. A projectile is accelerated in a handpiece and hits an applicator, which converts the mechanical energy into radial pressure waves. These propagate in a divergent manner and primarily affect superficial tissue layers. The energy decreases with increasing penetration depth, so that the therapeutic effect is generally limited to tissue at a depth of up to two centimetres. Due to the large-area distribution of energy, radial pressure wave therapy is particularly suitable for structures close to the surface such as tendon insertions, muscles and soft tissue. Typical areas of application are the treatment of muscle tension, trigger points, heel spurs or tennis elbow. It is also frequently used for achilles tendon irritation or myofascial pain syndromes.

Focused

In contrast to this is focused shock wave therapy, in which the waves are bundled and focused specifically on one point in the tissue. These waves are generated using piezoelectric, electromagnetic or electrohydraulic principles. This allows very high energy densities to be achieved, which penetrate the tissue up to a depth of ten centimetres at the focal point. Focused shock wave therapy enables precise treatment of deeper structures, such as deep-seated inflamed tendon attachments, calcium deposits on the shoulder or poorly healing bone fractures, so-called pseudarthroses. This method can also be used to treat complex pain syndromes that arise in deeper layers of tissue, such as the pelvic floor.

Another difference lies in the sensation of pain during treatment. Radial pressure wave therapy is generally perceived as painful. Focused shock wave therapy, on the other hand, can be pleasant, depending on the energy level selected and the depth of the tissue to be treated. Radial devices insert their pressure wave starting at the skin. This has many pain receptors which are irritated and therefore cannot - like the focused shock wave - have a precise effect on a pathological area. With focused devices, on the other hand, only the nociceptors/pain receptors that are "hypersensitive" due to pathology react to the mechanical pressure wave.

Both procedures are considered safe, provided they are carried out professionally. Side effects are usually minor and include temporary reddening of the skin, swelling or bruising.

Ultimately, the choice of the appropriate form of therapy depends on the individual diagnosis and the depth of the structure to be treated. While radial pressure wave therapy has its limitations for superficial complaints, focused shock wave therapy offers advantages for deeper or chronic conditions. A comprehensive consultation with experienced specialists is recommended in all cases.

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