Squeezi®

1. Gebruiksaanwijzing Nederlands
2. Instructions for use English
3. Gebrauchsanweisung German
4. Mode d'emploi French
5. Istruzioni per l'uso Italian
6. Instrucciones de uso Spanish
7. Podręcznik użytkownika Polish

1. NEDERLANDS

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2. ENGLISH

WARNING: Read carefully and completely the Squeezi User Manual before using it. This User Guide contains important safety and health information. Additionally, it is also helpful to obtain additional information about Creative Therapy's software with which to use the Squeezi.

Do not use unlicensed accessories that may cause injury to yourself and others as well as performance degrading results to the system.

All Creative Therapy logos are trademarked by Creative Therapy.

© 2022 Creative Therapy.

2.1. Safety Guidelines

Prior to use, check that Squeezi is in good condition. Squeezi's carry the CE mark and meet the following requirements:

93/42/EEC for Class 1 medical devices - (EU) 2017/745

PPE1, Directive 89/686/EEC: for activities of minimal risk - (EU) 2016/425

Important precautions when Squeezi is used by children, the elderly or other vulnerable populations. An adult must always communicate all safety and health information to all children who wish to use Squeezi. Failure or failure to do so may result in injury or damage to the system or other property. Observe Squeezi users to ensure proper and safe use of the Squeezi.

An adult should supervise a young child when first using the Squeezi to be sure they understand and will use it properly. To avoid injury and damage to the Squeezi or surrounding objects, it is advised to follow the following guidelines during use:

Be attentive to your balance so that your Squeezi does not hit anything.

If your hands would get sweaty while playing games, stop playing for a while and dry your hands as well as the Squeezi before continuing to play.

- Make sure no more than 1 person is using Squeezi at a time.
- Make sure there is enough space between yourself and other objects and people while playing a game. One and a half meters is advised.
- Avoid stepping and jumping on Squeezi. Do not throw Squeezi.
- Please do not use Squeezi in wet and humid environments.
- Do not use Squeezi in environments below 10°C or above 45°C.
- Do not use Squeezi when you are under the influence of alcohol or drugs.
- Clean Squeezi with a damp cloth when dirt sticks to the surface.
- Use Squeezi only indoors.

Playing with Squeezi may require energy or physical activity. Consult a doctor before starting to use the Squeezi if:

- You have heart, respiratory, back or other orthopedic problems;

- You have high blood pressure;
- You are prohibited from physical activity.

If you are using the Squeezi and you should experience fatigue, shortness of breath, chest pain, dizziness or discomfort due to pain then stop playing immediately and consult a doctor as soon as possible.

2.2 Indications

Squeezi is an interactive movement tool. It serves as a motivator and mobilizer in rehabilitation therapy, entertainment, active lifestyle, elderly care, care of persons with disabilities, patients, persons with cerebral palsy, neuromuscular disorders, neurological disorders, persons with post traumatic injury, fall prevention, dementia and other cases where a movement tool is desirable. Note: Squeezi is not certified as a medical measurement tool.

2.3. Purpose of use

Squeezi serves as a motivational movement tool. They are intended for children aged 6+ years to elderly people aged 100+ years. The intention is to use Squeezi in conjunction with a smartphone, tablet, computer or television. Squeezi should be used under the supervision of a therapist, or in children of an adult.

The most common movements with Squeezi are described below. Other movements may be possible depending on the type of game that will be played and the creative interpretation of the facilitator.

- Squeezing with full hand
- Squeezing with fingers individually
- Reorienting your upper body in relation to the game
- Moving your upper body according to the rhythm of the music
- Moving Squeezi from point A to point B
- Moving the arm or hand away from your body
- Moving the arm or hand toward your body
- Pronation or supination of the palm of the hand
- Shoulder, elbow or wrist rotations

2.4. Installation

To operate Squeezi, the installation of the Squeezi app is required. This can be found in the Play Store or App Store of the device you want to install it on. You can install the app on both Android and iOS. You can install the app on your smartphone, laptop, desktop computer or tablet. Make sure the Bluetooth reception of this device is on.

Before the first use, it is recommended that Squeezi is fully charged. You do this by turning Squeezi on, this causes the light next to the switch to turn on. Then you plug the

included cable into Squeezi as well as a power source for at least four hours. Through the app, you will be able to consult Squeezi's battery level.

Squeezi is best used immediately right in front of a television, computer or tablet screen. The ideal distance from the screen is from 1.5m to 4m, depending on the size of the screen.

Follow the guidelines below when using the Squeezi:

- Make sure there is enough space between yourself and other objects or people during that the game is being played.
- Use the Squeezi horizontally on a stable surface, avoid slippery surfaces so that the user cannot fall.

Important: Squeezi will not work until the correct software is installed on the computer or tablet.

2.5. Synchronizing the Squeezi with the Creative Therapy software

To ensure proper operation of the Squeezi, it must be linked to the Creative Therapy software. Follow the steps below to connect the Squeezi. This procedure should be followed each time Squeezi is started.

- Turn on the Squeezi via the slider.
- When it is on, you will see a light next to the slider
- Open the Creative Therapy Software called 'Squeezi'.
- Connect Squeezi and the software by clicking on 'connect
- Squeezi will now connect to your software
- Games or analyses can now be started

The basic operations using Squeezi are done with physical positions as already described above. Other modes of gameplay may be used depending on the game chosen. Review the instructional explanation before playing the game. Also review the health and safety information on the previous pages and the instructional explanation.

2.6. Lifetime

For optimal service life, the product must be handled carefully and cautiously. The flexible parts of Squeezi must be replaced with new ones after 2 years, since optimal flexibility can no longer be guaranteed at that time due to aging materials.

The product must be replaced based on one or more of the following criteria:

- Shape changes, e.g. due to intensive use or misuse
- Exceptional overloading
- Frequency of use
- Degree to which the product has been cared for
- Hygienic conditions of Squeezi

2.7. Squeezi cleaning and care

Squeezi can be cleaned with mild soap or mild disinfectant and a damp cloth or sponge. Avoid water and strong solvents. Strong solvents may contain substances that may damage Squeezi and thus adversely affect the performance of the product. Squeezi should not be washed. Do not iron Squeezi or expose it to direct heat sources such as stoves, heaters, radiators, direct sunlight, etc.

2.8. Accessories

Various options are available on Creative Therapy products at an additional cost. Please contact your local specialist or Squeezi manufacturer for more information.

2.9. Troubleshooting

If Squeezi refuses to function, or does not work as expected, follow the steps below before contacting help.

Squeezi is not responding:

- Make sure Squeezi is turned on
- Make sure Squeezi is charged
- Make sure Squeezi is properly connected to the software
- Avoid interference with other devices. If Bluetooth communication appears to be affected by other devices, please turn them off or move Squeezi to another room. Devices that may interfere with Bluetooth are: WIFI, microwave ovens, smartphones and other wireless devices.
- Review the synchronization roadmap and make sure all steps have been followed correctly. Once all steps have been followed, the following can be tried:
 - Turn off Squeezi. Wait at least 30 seconds, turn Squeezi back on and wait 1 minute.
 - Exit and restart the Creative Therapy software.
 - Check the FAQ page at www.creativetherapy.be. A similar problem may be described there.

If none of the above solutions work, you can always contact us at support@mycreativetherapy.com.

2.10. Warranty

This product is covered by a two-year warranty. The warranty period starts on the date of purchase (see receipt). During the warranty period you can submit complaints about a defective Squeezi to your dealer or to the specified service address. During this period, we will repair all defects resulting from defective material or operation at no charge. Defects resulting from careless handling, falls, collisions, water damage and the like are excluded from warranty. If Squeezi has been modified or altered and therefore is no longer in its original condition, the warranty will be void. If Squeezi was repaired or

exchanged, the old warranty period will be maintained; the warranty period of 2 years, starting from the date of purchase, remains unchanged. The warranty expires in case of non-compliance, such as use other than the intended use, and repairs and modifications of Squeezi or parts by unauthorized persons.

Thank you for your interest in a Creative Therapy product.

If you have any questions or need assistance, please feel free to contact us or our dealers and representatives.

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