


Flip and go!




RIWObutton

Implant for femoral fixation
of ligament grafts

RIWObutton

Implant for femoral fixation of ligament grafts

- Especially for the fixation of ligament grafts for crucial ligament reconstruction
- The thread loop is supplied in lengths from 15 – 35 mm in 5 mm increments. This permits optimum adjustment of the implant to the anatomical requirements of the specific patient and simplifies preparation.
- No change in length of the loop
- Fast transplant preparation due to preconfigured loop on titanium button
- Broad contact surface between transplant and loop
- Low cost



	Loop length	Type
RIWObutton sterile, for tendon fixation, 2 threads (green/white), Ø 0.7 mm	15 mm	491801015
	20 mm	491801020
	25 mm	491801025
	30 mm	491801030
	35 mm	491801035