



Fluid Control 2203

Irrigation pump for Arthroscopy

- Non-invasive pressure measurement
- Optimum operating convenience
- CAN-Bus interface
- High efficiency
- Instrument recognition
- Improved flow rate



Fluid Control 2203

Irrigation pump for Arthroscopy

Fluid Control 2203

Irrigation pump for Arthroscopy

- Non-invasive pressure measurement of fluids
- Clear display provides optimum operating status
- CAN-Bus interface for voice control
- Maximum convenience provided by separate daily and single use patient tubing sets
- Instrument recognition
- Improved flow rate in narrow-lumen instruments

Fluid Control 2203

and mains cable,
length of 3 m (2440.03).....2203.001

Alternative

with Can-Bus interface.....2203.002



Patient-specific set

with Safe Lock, disposable,
PVC, sterile, consisting of
1 x daily tubing sets 4170.229 and
10 x patient-specific tubing sets ..4170.228

Daily tubing sets

with Safe Lock, disposable,
PVC, sterile, 10 pieces.....4170.229

Foot switch

pneumatic2030.231

Remote control cable

autoclavable.....2202.961
non autoclavable.....2202.991

Technical data

Power supply100-240 V, 50/60 Hz
Power consumptionmax. 50 W
Safety class according to VDE
0750 / IEC 6011
ClassificationBF
Pressure range
(infinitely adjustable)0-150 mm Hg
Dimensions
(W x H x D).....360 x 150 x 350
Weight.....8.4 kg
Can-Bus interfaceonly for 2203.002
Software update for instrument recognition

Accessories

Patient-specific set

with piercing spike, disposable,
PVC, sterile, consisting of
1 x daily tubing sets 4170.226 and
10 x patient-specific tubing sets ..4170.227

Daily tubing sets

with piercing spike, disposable,
PVC, sterile, 10 pieces4170.226