

Access. Diagnose. Treat.



Viper™

Flexible Ureteroscope

Innovations in Endourology from Richard Wolf

Viper™



The Trusted, Long-Standing Solution


The Viper has an award-winning ergonomic design that boasts a light-weight, easy-to-hold body that fits almost every hand size and grip.

What Are the Surgical Benefits?

The Viper's small 6 Fr. distal tip self-dilates to 8.8 Fr., allowing easy introduction, increased durability, and stability.

Utility and Flexibility

The working channel accommodates numerous 3 Fr. accessory devices, including graspers, biopsy forceps, and electrodes.



The Go-To Solution for Flexible Ureteroscopy

Viper™ Flexible Ureteroscope



Visualization

- 270° up/down distal deflection allows complete visualization of the kidney
- Continuous cyclonic flow clears blood and debris for a crystal-clear image
- State-of-the-art fiber optic bundle provides superior images
- 680 mm length is ideal for access to all areas of the upper urinary tract

Versatility

- Large 3.6 Fr. working channel allows for a variety of accessory instruments
- Atraumatic 6 Fr. distal tip dilates to 8.8 Fr. for ease of entry
- Award-winning ergonomic handle reduces fatigue and increases comfort
- Built-in 3 mm luer bridge facilitates easy connections
- Compatible with numerous reprocessing options, including: EO, STERIS System 1E®, and STERRAD®



Dependability

- Proven durability allows for long-term use of product
- Flexibility permits better access to hard-to-reach areas within the upper urinary tract
- Self-dilating distal tip combines small size without sacrificing utility

Ordering Information

Viper Flexible Ureteroscope	
Description	Part Number
Viper Flexible Ureteroscope Kit, Ergonomic Deflection	7325.0761
Viper Flexible Ureteroscope Kit, Logical Deflection	7325.0711
2.5 mm Armored Light Cable	8061.256
Sealing Adapter, 0-6 Fr. Capacity, Sterile, 5 per box	20.SA65
3 Fr. Flexible Biopsy Forceps, 920 mm WL	829.601
3 Fr. Flexible Grasping Forceps, 920 mm WL	828.651
3 Fr. Coagulating Electrode, 920 mm WL	8741.80
Tray for Flexible Fiberscopes	85850302



Visit our website for more information about the Richard Wolf endourology solutions, or contact your local representative to schedule an evaluation!



Technical specifications subject to change without notice.

971-02.01-0113USA

Financing and service agreements are available. Please contact your local Richard Wolf representative for details.
 353 Corporate Woods Parkway • Vernon Hills, Illinois 60061-3110 • Phone: (800) 323-WOLF (9653) • Fax: (847) 913-1488
 Email: sales@richardwolfusa.com • Website: www.richardwolfusa.com