


Dual uniqueness
meets digital vision



spirit of excellence



Watch our trailer.

BOA
VISION

&


COBRA
VISION

Flexible single and dual channel Sensor Uretero Renoscopes



Flexible single channel Sensor Uretero Renoscope

- Fits in access sheaths from 9.5 Fr.
- LED illumination and distal sensor
- Highest torsional stiffness







Flexible dual channel Sensor Uretero Renoscope



spirit of excellence

- Slim laser channel with laser shifter
- Large working and irrigation channel
- LED illumination and distal sensor
- Highest torsional stiffness



R. Wolf Sensor URS & ENDOCAM® Logic HD

Dual uniqueness meets digital vision



spirit of excellence

"BOA VISION" – Flexible single channel Sensor Uretero Renoscope

distal sensor, integrated LED illumination, 6 / 8.9 Fr., DEF 270° up/
down, passive DEF, WL 680 mm,

control lever to distal – downward DEF 7355071

control lever to distal – upward DEF 7355076

"COBRA VISION" – Flexible dual channel Sensor Uretero Renoscope

distal sensor, integrated LED illumination, 6 / 9.9 Fr., DEF 270° up/
down, passive DEF, WL 680 mm, with special laser shifter,

control lever to distal – downward DEF 7356071

control lever to distal – upward DEF 7356076



ENDOCAM® LOGIC HD Controller Set

comprising: ENDOCAM® LOGIC HD Controller, manual remote
control, USB flash drive 8 GB, CAN BUS Terminator 55251071

ENDOCAM® LOGIC HD Lite Controller Set

comprising: ENDOCAM® LOGIC HD Lite Controller,
manual remote control, USB flash drive 8 GB 55252011

www.richard-wolf.com