

USB 300



USB 300

The low cost HD recording solution

- Records videos and images in true 1080p HD
- Intuitive controls for simplified use
- Data storage on USB flash drive or on internal hard drive with network accessibility

USB 300

Intuitive HD recording at an affordable price

External Storage

Directly export HD images and videos on a USB flash drive



Internal Storage

HD images and videos will be saved on 320 GB internal hard drive and can be accessed through network.

Complete your HD experience with Richard Wolf's fully autoclavable 5550 HD Camera.



Technical Specification:

LCD Monitor	2 1/2" full color (for video confirmation only, not for diagnosis)
Buttons	Sealed membrane, fluid resistant
Video Input	DVI (RGB via adapter), HD SDI, S-Video, Composite
Video Output	DVI (RGB, VGA via adapter)
Video Formats	H.264 / MPEG4 (1080, 720, NTSC, PAL)
Image Formats	JPG, BMP
Networking	RJ45 10/100/1000 Ethernet
Internal Storage	Internal hard drive, 320GB (for use with networking systems only)
External Media	USB flash drives or USB hard drive
USB Support	USB2.0, high speed, 1 front connector, 2 rear connectors
Inputs	Standard 3.5 mm mini jacks (one for photos, one for videos)
Dimensions	240 x 211 x 63 mm
Weight	1.8 kg / 3 lb 15 oz
Power Supply	AC 100-240V 50/60Hz 25W 3-prong AC jack
Certifications / Compliance	CE 93/42/EEC UL 60601-1, IEC 60601-1, IEC 60601-1-2 CSA C22.2 No.601.1, E 306417 FCC Part 15, FDA Class 1, RoHS

Ordering Information:

USB 300 HD Recorder.....	5654.007
includes: 4GB Flashdrive	5654.0044
Keyboard.....	41065.5654
Video and Power Cables	

For more information contact your local Richard Wolf representative or call Customer Service at (800) 323-WOLF (9653)

353 Corporate Woods Parkway • Vernon Hills, Illinois 60061-3110 • Phone: (800) 323-WOLF (9653) • Fax: (847) 913-1488
Email: sales&marketing@richardwolfusa.com • Website: www.richardwolfusa.com

