

Baseline[®] volumetric edema gauges

- objectively measure body part volume by using the fluid displacement method
- body part is placed in a water-filled tank; the water overflow is collected and measured
- amount of edema can be measured and quantified



each set includes

- measurement tank
- overflow beaker
- graduated cylinder


upper extremity (LxWxH)

12-3500	hand set 5" x 3" x 9"	200.00
12-3501	arm set 7" x 7" x 30"	300.00
12-3504	forearm set 6" x 6" x 24"	250.00

lower extremity (LxWxH)

12-3502	foot set 13" x 5" x 9"	225.00
12-3505	partial leg set 13" x 6" x 24"	350.00
12-3508	large partial leg set 13" x 8" x 16"	365.00

Baseline[®] measuring tape



- reads in cm and inches
- locks in place until reset button is pressed
- 3 lengths available

12-1210	5' tape (each)	5.00
12-1210-25	25 each	105.00
12-1211	6' tape (each)	8.00
12-1211-25	25 each	170.00
12-1212	10' tape (each)	12.50
12-1212-25	25 each	265.00

Gulick measurement tape

- spring exerts constant force on tape for accurate and repeatable results
- push-button retractor
- reads in cm and inches



12-1201	5' (each)	22.50
12-1201-25	25 each	477.50
12-1203	6' (each)	27.50
12-1203-25	25 each	585.00
12-1204	10' (each)	32.50
12-1204-25	25 each	690.00

Collyer breadth caliper

- breadth caliper can be used to measure chest, pelvis, shoulder and more
- 15" scale is marked in both cm and inches

12-1231	Collyer caliper	30.00
---------	-----------------	-------

Martin breadth caliper

- for accurate breadth readings for chest, pelvis, shoulder, and more
- 20" scale is marked in both cm and inches

12-1230	Martin caliper	35.00
---------	----------------	-------

stethoscopes

nurse's

- anodized aluminum chest piece, plastic ear tips, 22" tube



12-2200	nurse's (each)	7.50
12-2200-25	25 each	157.50

dual head

- turn to convert to non-child ring; external spring binaural, 22" tube



12-2210	dual-head (each)	10.00
12-2210-25	25 each	235.00
12-2211	stainless steel, latex free	52.50

Sprague-Rappaport

- rotating drum stem and external spring binaural: with 3 bell chest pieces and 2 diaphragms, 22" tube



12-2220	Sprague-Rappaport	25.00
---------	-------------------	-------

Baseline[®] circumference tape

- the circumference of a body part measure in both cm and inches
- measure by yourself: simply make loop, position loosely around waist, wrist or area being measured, and press button
- measuring tape will stay in place securely while you have two free hands to make adjustments and record measurements



Baseline[®] finger circumference gauge

- wrap webbing around finger or digit and read circumference in both cm and inches



12-1205	circumference tape (each)	10.00	12-1222	gauge in/cm	12.50
12-1205-25	25 each	212.50	12-1222-25	25 each	295.00

evaluation